<b>CACIPLIQ20® Spray 7.5 RGTA®</b> Cutaneous chronic wound healing solution Matrix therapy for regeneration	EN
---	----

#### INSTRUCTIONS FOR USE

##### ( Composition )

Sodium salt of a 1-6 poly(carboxymethylglucose sulphate) 0.1 mg/mL, sodium chloride 9 mg/mL, sodium hydroxide, purified water.

##### ( General Information )

CACIPLIQ20® is an RGTA® (ReGeneraTing Agents) extracellular matrix mimetic agent, which promotes natural skin regeneration. It is supplied in a 7.5 mL spray bottle of CACIPLIQ20® solution. Each spray corresponds to a dose of 14 µg of RGTA®. CACIPLIQ20® is a solution intended to improve the healing of chronic cutaneous ulcers by promoting the restoration of the extracellular matrix, with the consequent reduction in lesion-induced pain.

##### ( Indication )

CACIPLIQ20® is indicated for chronic cutaneous ulcers caused by pressure (bedsores), peripheral arterial disease (e.g. stage IV Leriche and Fontaine classification) or diabetes (including amputation), showing no tendency to heal after 6 months of standard care, or still unhealed after 12 months.

##### ( User instructions )

Normal hygiene precautions should be taken when treating a wound with CACIPLIQ20® (e.g. handwashing, appropriate protective clothing, examination table sheet...). In case of a sharp debridement, follow current best practice procedures to ensure sterility and prevent cross contamination. After cleaning, disinfecting, rinsing with sterile physiological saline serum and drying the wound, spray the surface and edges of the wound from 5 cm away. Wait 5 minutes and then dress the wound with the usual gel dressing and a non-adherent bandage. CACIPLIQ20® is intended to be applied twice a week (every 3 to 4 days), ideally when the dressing is changed, with only one spray per affected area (7-9 cm<sup>2</sup>). Results may vary from patient to patient, your physician may adapt the duration of treatment accordingly. Results can generally be seen within 3 months of use.

##### ( Contraindications )

Do not use on patients with a history of allergy or anaphylactic reactions to heparin or other heparinoids. Individuals sensitive to heparinoids may have a local reaction, (redness, local swelling, pain, etc.), which should not affect the healing process. If symptoms persist, or worsen, discontinue use and seek medical advice. Due to lack of clinical data, CACIPLIQ20® is not indicated for use by pregnant or breastfeeding patients; nor patients under 18 years old. CACIPLIQ20® is not indicated for patients eligible for revascularization surgery. CACIPLIQ20® is not indicated for patients using a non-removable discharge boot to treat diabetic foot ulcers. Efficacy of CACIPLIQ20® may be reduced on necrotic tissue or purulent wounds. It is not recommended to use CACIPLIQ20® in combination with locally applied aminoglycoside antibiotics (gentamicin, neomycin, etc.), solutions containing silver salts, copper, gold or iodine, or polycationic-based products such as quaternary ammonium disinfectants (e.g. benzalkonium chloride).

##### ( Warnings )

Respect the dosage prescribed. If CACIPLIQ20® is applied too frequently, it may reduce the therapeutic effect. Do not inject. Do not swallow. Keep out of the reach and sight of children. Do not re-sterilize. The product is sterile for up to 6 months after first use. However, if the spray head comes into contact with the wound, or is otherwise contaminated, disinfect external surface of the pump before the next use (e.g. with an alcohol wipe).

##### ( Precautions )

Do not use if the packaging is damaged. In the event of unintended contact with the eye, rinse with lukewarm running water. To reduce the risk of cross contamination, it is not recommended to use a single bottle of CACIPLIQ20® on multiple patients.

##### ( Users )

CACIPLIQ20® should preferably be applied to a wound cleaned by sharp debridement by a health care professional. CACIPLIQ20® should be used under the supervision of a health care professional.

##### ( Storage )

Store in a dry place between 2°C and 25°C. Brief and occasional exposure to higher temperatures (less than 40°C) does not harm the integrity of the product. Do not use after the expiry date indicated on the pack. Do not use if the bottle or spray head is damaged and not sterile. Use within 6 months of the first spray.

##### ( Recycling )

After use, dispose of: the spray in a glass recycling bin and the cardboard packaging in a cardboard recycling bin.

Please refer to this link to find the SSCP (Summary of Safety and Clinical performance) on EUDAMED database: <https://ec.europa.eu/tools/eudamed>

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